

# S1 TITAN/TRACER 5/CTX

- **EasyCal— A Comic Strip**


**V.2.5.61**

# A COMIC STRIP





---

## Contents

1	Setting up the Workspace.....	4
2.	A New Calibration. ....	5
2.1	Select the elements. ....	6
2.2	Adding Compounds.....	8
2.3	Enter Standard Sample Names and Concentrations.....	9
3	Measuring the Assays. ....	12
3.1	Creating the Run Order.....	12
3.2	Transferring the Run Order to the Instrument.....	13
3.3	Running the Run Order.....	14
3.4	Transferring the Data to the PC.....	15
4.	Setting Envelop background subtraction.....	19
5.	Calibration.....	20
5.1	Importing the data into EasyCal.....	20
6.	Modules.....	34
7.	Saving the calibration.....	35
7.1	The appinstall package.....	35
7.2	Installing the appinstall package.....	36
7.3	Testing the calibration.....	37




# 1 Setting up the Workspace.


 EasyCal Launcher			<h2 style="text-align: center;">S2Setup</h2>  <p>See Appendix A for more information about "Project Base Folder"</p>
Double Click			


## 2. A New Calibration.

<p>Note</p>	<p>The internal standard selected when creating a new calibration will be used as the default for each element as they are added to the calibration. The active internal standard can be changed on an element by element basis at a later point in the process. The default selected at this point cannot be changed.</p>




## 2.1 Select the elements.





## 2.2 Adding Compounds.





## 2.3 Enter Standard Sample Names and Concentrations.

The first screenshot shows the 'Material' field in the 'Material group' section. A red box highlights the text 'SS401' in the 'Material' dropdown menu, with a green arrow pointing to it from a text box that says 'Enter the name of Standard Samples'.

The second screenshot shows a table of standard samples. A green box highlights the table, with a text box that says 'Enter all Standard Sample Names'. Below the table, a green box highlights the 'Custom Units' field, with a text box that says 'Units: PPM, % or Custom'.

	Sum(%)	V(%)	Cr(%)	Mn(%)	Fe(%)	Ni(%)	Cu(%)	Mo(%)
SS401	0.00	0.496	T	T	T	T	T	T
SS402	0.00	0.194	T	T	T	T	T	T
SS403	0.00		T	T	T	T	T	T
SS404	0.00		T	T	T	T	T	T
SS405	0.00		T	T	T	T	T	T
SS406	0.00		T	T	T	T	T	T
SS407	0.00		T	T	T	T	T	T
SS408	0.00		T	T	T	T	T	T
SS409	0.00		T	T	T	T	T	T

Material group: **Stainless**

Limits: 1.0 105.0 %

After validation, go to: Bottom Right

	Sum(%)	V(%)	Cr(%)	Mn(%)	Fe(%)	Ni(%)	Cu(%)	Mo(%)
SS401	0.50	0.496	T	T	T	T	T	T
SS402	0.19	0.194	T	T	T	T	T	T
SS403	0.00	T	T	T	T	T	T	T
SS404	0.00	T	T	T	T	T	T	T
SS405	0.00	T	T	T	T	T	T	T
SS406	0.00	T	T	T	T	T	T	T
SS407	0.00	T	T	T	T	T	T	T
SS408	0.00	T	T	T	T	T	T	T
SS409	0.00	T	T	T	T	T	T	T

Default is Trace

Material group: **Stainless**

Limits: 1.0 105.0 %

After validation, go to: Bottom Right

	Sum(%)	V(%)	Cr(%)	Mn(%)	Fe(%)	Ni(%)	Cu(%)	Mo(%)
SS401	2.45	0.496	0.138	1.197	T	0.019	0.101	0.495
SS402	2.32	0.194	0.652	0.228	T	0.808	0.302	0.14
SS403	3.01	0.341	0.463	1.677	T	0.223	0.221	0.088
SS404	2.54	0.107	0.774	0.532	T	0.393	0.427	0.307
SS405	1.67	0.411	0.206	0.903	T	0.102	0.022	0.025
SS406	5.35	0.01	2.001	0.447	T	1.62	0.289	0.98
SS407	5.17	0.19	3.03	0.195	T	0.527	0.397	0.83
SS408	5.66	0.067	0.111	0.557	T	4.13	0.694	0.098
SS409	5.71	0.008	1.318	0.559	T	3.02	0.205	0.599

-OR-

Paste with Titles

Material group: **Stainless**

Limits: 1.0 105.0 %

After validation, go to: Bottom Right

	Sum(%)	Ti(%)	V(%)	Cr(%)	Fe(%)	Ni(%)	Cu(%)	Mo(%)
SS401	1.25	T	0.496	0.138	T	0.019	0.101	0.495
SS402	2.10	T	0.194	0.652	T	0.808	0.302	0.14
SS403	1.34	T	0.341	0.463	T	0.223	0.221	0.088
SS404	2.01	T	0.107	0.774	T	0.393	0.427	0.307
SS405	0.77	T	0.411	0.206	T	0.102	0.022	0.025
SS406	4.90	T	0.01	2.001	T	1.62	0.289	0.98
SS407	4.97	T	0.19	3.03	T	0.527	0.397	0.83
SS408	5.10	T	0.067	0.111	T	4.13	0.694	0.098
SS409	5.15	T	0.008	1.318	T	3.02	0.205	0.599

Material group: **Stainless**

Limits: 1.0 105.0 %

After validation, go to: Bottom Right

	Sum(%)	V(%)	Cr(%)	Mn(%)	Fe(%)	Ni(%)	Cu(%)	Mo(%)
SS401	2.45	0.496	0.138	1.197	T	0.019	0.101	0.495
SS402	2.32	0.194	0.652	0.228	T	0.808	0.302	0.14
SS403	3.01	0.341	0.463	1.677	T	0.223	0.221	0.088
SS404	2.54	0.107	0.774	0.532	T	0.393	0.427	0.307
SS405	1.67	0.411	0.206	0.903	T	0.102	0.022	0.025
SS406	5.35	0.01	2.001	0.447	T	1.62	0.289	0.98
SS407	5.17	0.19	3.03	0.195	T	0.527	0.397	0.83
SS408	5.66	0.067	0.111	0.557	T	4.13	0.694	0.098
SS409	5.71	0.008	1.318	0.559	T	3.02	0.205	0.599

- Copy
- Copy All
- Paste
- Trace
- Inaccurate
- Column Balance
- Line Balance
- Unit
- Format
- Column Decimals

Material group: Inconel

Limits: 1.0 105.0 %

After validation, go to: Bottom Right

	Sum(%)	V(%)	Cr(%)	Mn(%)	Fe(%)	Ni(%)	Cu(%)	Mo(%)
SS401	100.00	0.496	0.138	1.197	97.554	0.019	0.101	0.495
SS402	100.00	0.194	0.652	0.228	97.6			
SS403	100.00	0.341	0.463	1.677	96.9			
SS404	100.00	0.107	0.774	0.532	97			
SS405	100.00	0.411	0.206	0.903	98.3			
SS406	100.00	0.01	2.001	0.447	94.6			
SS407	100.00	0.19	3.03	0.195	94.8			
SS408	100.00	0.067	0.111	0.557	94.3			
SS409	100.00	0.008	1.318	0.559	94.2			

- Copy
- Copy All
- Paste
- Trace
- Inaccurate
- Column Balance
- Line Balance
- Unit >
- Format >
- Column Decimals >

Material group: Inconel

Limits: 1.0 105.0 %

After validation, go to: Bottom Right



	Sum(%)	V(%)	Cr(%)	Mn(%)	Fe(%)	Ni(%)	Cu(%)	Mo(%)
SS401	100.00	0.496	0.138	1.197	97.554	0.019	0.101	0.495
SS402	100.00	0.194	0.652	0.228	97.676	0.808	0.302	0.14
SS403	100.00	0.341	0.463	1.677	96.987	0.223	0.221	0.088
SS404	100.00	0.107	0.774	0.532	97.46	0.393	0.427	0.307
SS405	100.00	0.411	0.206	0.903	98.334	0.102	0.022	0.025
SS406	100.00	0.01	2.001	0.447	94.653	1.62	0.289	0.98
SS407	100.00	0.19	3.03	0.195	94.834	0.527	0.397	0.83
SS408	100.00	0.067	0.111	0.557	94.343	4.13	0.694	0.098
SS409	100.00	0.008	1.318	0.559	94.291	3.02	0.205	0.599



ApplicationWizard

The definition of the material has been changed at this step; do you want to save the changes?  
Note: the standard materials will be saved in their modified state as well.


### 3 Measuring the Assays.

#### 3.1 Creating the Run Order.






### 3.2 Transferring the Run Order to the Instrument.





Drag csv file to USB stick into RunOrders folder




Move USB stick to the instrument

### 3.3 Running the Run Order.





### 3.4 Transferring the Data to the PC.

Move USB stick to computer





Drag the 7z file from the USB to the computer.















## 4. Setting Envelop background subtraction

	<p>The Envelope background option found on the Calibration Menu can be used to configure how element peak intensities are calculated. Use of Net intensity, the background under the peak is subtracted, is the default. This behavior can be changed by selecting options from the Envelope background parameters form.</p>
	<p>Uncheck the “Use net intensities in empirical model” check box to use Raw peak intensities.</p> <p>* Checking the “Use those settings as default” box and clicking the “OK” button will cause the entered values to be used as the default for all calibrations created subsequent to this change.</p> <p>** The other values found on this form should only be changed if the user is an expert with this tool.</p> <p>*** CalToolkitEC uses 2.0 and 1.6 as well.</p>

## 5. Calibration.


### 5.1 Importing the data into EasyCal.



	<p>Select the element with the largest Z value. Click on the two left/right arrow buttons to scroll through the elements. Or click on the button with the element, which will bring up the periodic chart.</p>


**Internal Standards**

Note: Each element in the calibration can use a different internal standard.  
A single range can be defined for the User ROI option.




Internal Standard:  
**None**

This selection uses the assay acquisition length to normalize the peak intensity to counts per second.




Internal Standard:  
**Sum of Intensities**

This selection uses the accumulated peak intensities of all other elements in the calibration as the normalizing divisor.





**Internal Standard:  
Compton (Rh)**

This selection uses the Compton scattering region of the x-ray source target element as the normalizing divisor. Generally this will be the Compton region of the Rh element.




**Internal Standard:  
User ROI**


This selection is similar to Compton (Rh) with the exception that the user sets the ROI that's used as the normalizer. The region is defined using the Min and Max fields. The values are entered in keV of energy.








Select the columns to be displayed in the <Standards> panel




Select the next element (Cu) and click on the <Standards> tab. Note the "Std Dev."

Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation/	Relative Deviation	Gross Intensity
9	SS409_Solids	0.2050	0.2703	0.0653	32	9.005
4	SS404_Solids	0.4270	0.3980	-0.0290	-9.1	13.099
3	SS403_Solids	0.2210	0.1937	-0.0273	-12	6.822
2	SS402_Solids	0.3020	0.2830	-0.0190	-6.3	9.759
6	SS406_Solids	0.2890	0.3062	0.0172	6.0	10.031
7	SS407_Solids	0.3970	0.3843	-0.0127	-3.2	12.152
8	SS408_Solids	0.6940	0.7052	0.0112	1.6	23.904
5	SS405_Solids	0.0220	0.0302	0.0082	37	1.307
1	SS401_Solids	0.1010	0.0971	-0.0039	-3.9	3.546



Sort the <Absolute Deviation> column and note that SS409 is the largest.



Select the <Line Over.> tab. Click on <Insert New...> or right-click on blank space to insert new.



Select the best overlap suggestion (Ni KB1 in this case).



Note the line overlap is added to the table. To delete right click on the entry. Select the <Standards> tab.

Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation/	Relative Deviation	Gross Intensity
7	SS407_Solids	0.3970	0.4165	0.0195	4.9	12.152
3	SS403_Solids	0.2210	0.2067	-0.0143	-6.5	6.822
8	SS408_Solids	0.6940	0.6861	-0.0079	-1.1	23.904
2	SS402_Solids	0.3020	0.2941	-0.0079	-2.6	9.759
9	SS409_Solids	0.2050	0.2121	0.0071	3.5	9.005
6	SS406_Solids	0.2630	0.2541	-0.0089	-1.8	10.881
4	SS404_Solids	0.4270	0.4261	-0.0009	-0.20	13.099
5	SS405_Solids	0.0220	0.0214	-0.0006	-2.6	1.307
1	SS401_Solids	0.1010	0.1009	-0.0001	-0.12	3.546

Note the <Absolute Deviation> is now 0.0071. Less than 1/3 of the “Std Dev.”

Calibration Toolbox

**Ni**

Ni K $\alpha$ 1 Tutorial

StdDev: 0.062 %

Redo Undo

Offset Off On

Quadratic Off On


Calibration Plot for Ni

Standards


Next element (Ni). Note the “StdDev.” and select the <Standards> tab.

Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation/	Relative Deviation	Gross Intensity
8	SS408_Solids	4.1300	4.0296	-0.1004	-2.4	91.423
6	SS406_Solids	1.6200	1.7015	0.0815	5.0	36.707
9	SS409_Solids	3.0200	3.0935	0.0735	2.4	67.608
7	SS407_Solids	0.5270	0.5914	0.0644	12	12.386
1	SS401_Solids	0.0190	0.0607	0.0417	219	1.375
5	SS405_Solids	0.1020	0.1386	0.0366	36	3.139
3	SS403_Solids	0.2230	0.2494	0.0264	12	5.691
2	SS402_Solids	0.8080	0.7834	-0.0246	-3.0	17.747
4	SS404_Solids	0.3930	0.4066	0.0136	3.5	9.095

Note the <Absolute Deviation> is slightly larger than the “Std Dev.” However, when we go to the <Line Over.> tab there are no overlap suggestions.



The next element is Fe. Note that in the figure we clicked <Fix to current point s>to adjust the x-axis. Fe is regarded as the matrix. Because Fe is calculated as the balance, while ignoring minor elements such as P, Si, N, it is inaccurate. We do not calibrate Fe.



The next element is Mn. The “Std Dev.” is 0.33 %. The points are scattered.

Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation	Relative Deviation	Gross Intensity
7	SS407_Solids	0.1950	0.7088	0.5138	264	57.878
6	SS406_Solids	0.4470	0.7252	0.2782	62	60.272
8	SS408_Solids	0.5570	0.3253	-0.2317	-42	46.456
5	SS405_Solids	0.9030	0.6913	-0.2117	-23	61.834
1	SS401_Solids	1.1970	0.9917	-0.2053	-17	74.589
9	SS409_Solids	0.5590	0.6827	0.1237	22	59.325
3	SS403_Solids	1.6770	1.5536	-0.1234	-7.4	99.004
2	SS402_Solids	0.2280	0.1264	-0.1016	-45	37.977
4	SS404_Solids	0.5320	0.4900	-0.0420	-7.9	52.672

The <Absolute Deviation> is 0.5138 for SS407.

Overlap correction for Mn KA1/Tutorial (5.899 keV)

Overlapping line: **Cr KB1 (5.947 keV, 14%)** OK

Correct by: **Measured** **Cr KA1/Tutorial** Cancel

Factor:   Given  Adjustable

Search parameters


Min concentration:   Escape peaks

Distance (eV):   Sum peaks

In the <Line Over.> tab insert the overlap as shown. It is inserted in the table.

Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation	Relative Deviation	Gross Intensity
9	SS409_Solids	0.5590	0.5918	0.0328	5.9	59.325
5	SS405_Solids	0.9030	0.8844	-0.0186	-2.1	61.834
7	SS407_Solids	0.1950	0.1768	-0.0182	-9.4	57.878
2	SS402_Solids	0.2280	0.2135	-0.0121	-5.3	37.977
6	SS406_Solids	0.4470	0.4569	0.0099	2.2	60.272
4	SS404_Solids	0.5320	0.5397	0.0077	1.5	52.672
3	SS403_Solids	1.6770	1.6716	-0.0054	-0.32	99.004
8	SS408_Solids	0.5570	0.5593	0.0023	0.41	46.456
1	SS401_Solids	1.1970	1.1986	0.0016	0.13	74.589


The <Absolute Deviation> for SS407 is now down to 0.0182% and the largest one is still only 0.0328 %. Quite an improvement.



Next is Cr. Adding an offset improves the “Std Dev.” But let’s look at overlaps.

Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation	Relative Deviation	Gross Intensity
8	SS408_Solids	0.1110	0.0811	-0.0299	-27	21.85
5	SS405_Solids	0.2060	0.2187	0.0127	6.2	31.02
4	SS404_Solids	0.7740	0.7852	0.0112	1.5	68.08
1	SS401_Solids	0.1380	0.1486	0.0106	7.7	26.35
7	SS407_Solids	3.0300	3.0251	-0.0049	-0.16	202.39
6	SS406_Solids	2.0010	2.0032	0.0022	0.11	143.31
3	SS403_Solids	0.4630	0.4610	-0.0020	-0.43	47.59
2	SS402_Solids	0.6520	0.6513	-0.0007	-0.11	59.99
9	SS409_Solids	1.3180	1.3187	0.0007	0.051	100.97


The largest “Absolute Deviation” is 0.0299.



Select the V KB1 overlap, and click OK.

Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation	Relative Deviation	Gross Intensity
8	SS408_Solids	0.1110	0.0894	-0.0216	-19	21.85
4	SS404_Solids	0.7740	0.7898	0.0158	2.0	68.08
7	SS407_Solids	3.0300	3.0192	-0.0108	-0.36	202.39
9	SS409_Solids	1.3180	1.3251	0.0071	0.54	100.97
5	SS405_Solids	0.2060	0.2125	0.0065	3.2	31.02
6	SS406_Solids	2.0010	2.0075	0.0065	0.32	143.31
3	SS403_Solids	0.4630	0.4576	-0.0054	-1.2	47.59
1	SS401_Solids	0.1380	0.1390	0.0010	0.74	26.35
2	SS402_Solids	0.6520	0.6529	0.0009	0.13	59.99

It only improve the <Absolute Deviation" from 0.0299 to 0.0200.




Next select the Fe Escape peak. And click OK.

Number	Standard Name	Chemical Concentration	XRF Concentration	Absolute Deviation	Relative Deviation	Gross Intensity
8	SS408_Solids	0.1110	0.0947	-0.0163	-15	21.85
9	SS409_Solids	1.3180	1.3293	0.0113	0.86	100.97
7	SS407_Solids	3.0300	3.0190	-0.0110	-0.36	202.39
6	SS406_Solids	2.0010	2.0098	0.0088	0.44	143.31
4	SS404_Solids	0.7740	0.7827	0.0087	1.1	68.08
2	SS402_Solids	0.6520	0.6444	-0.0076	-1.2	59.99
1	SS401_Solids	0.1380	0.1443	0.0063	4.6	26.35
3	SS403_Solids	0.4630	0.4598	-0.0032	-0.70	47.59
5	SS405_Solids	0.2060	0.2090	0.0030	1.4	31.02

This again improves the deviations.

Overlap	Energy	Delta En.	Remark	Type	Base for calculation	Adjustable	Coefficient
Fe KB1	5.319	0.096	Escape	Measured	Fe KA1/Tracer	Adjustable	5.886e-005
V KB1	5.428	0.013	Measured	V KA1/Tracer	V KA1/Tracer	Adjustable	-0.0007717

We therefore remove the V KB1. Because there is so much more Fe, its affect outweighs the V affect.



Next is V. It is already good.

Tutorial.devm

ApplicationWizard

Save changes to Tutorial.devm?

Yes No Cancel

Save As


Save in: Calibrations

File name: Tutorial.devm


Save as type: ApplicationWizardDraft Files (\*.devm)

Save Cancel




Select the left most tab, <Tutorial.devm>. Then click . Click <Yes> to save the results.

## 6. Modules.



For Example:




Enter a name to show with results

Select a function

Enter the arguments


OR



Enter a name to show with results



Select a condition

Enter the arguments




## 7. Saving the calibration.


### 7.1 The appinstall package.



Check the <Chemistry Rules>




Select the <Info/Save> tab, then File -> Save to save the calibration installation package.




Note the folder as set by EasyCal where the appinstall package is stored.


## 7.2 Installing the appinstall package.



Open S1Sync, connect to the instrument. Open File Explorer and find the folder from 6.1. Drag and drop the appinstall file on the package folder.



### 7.3 Testing the calibration.



Login to the instrument and click <APPLICATION>	Select the new calibration and click <OK>	Note that the new name is displayed in the <APPLICATION> button	Take an assay and note the concentrations as well as the Sum and Ni/Cr ratio are shown
---	---	---	--

# The End